

Review Article





Sepsis-associated Acute Kidney Injury: Epidemiology, Pathophysiology, and Management Strategies in Critical Care

Ningyuan Xu and Linlin Zhang*

Department of Critical Care Medicine, Beijing Tiantan Hospital, Capital Medical University, Beijing, China

*Correspondence to: Linlin Zhang, Department of Critical Care Medicine, Beijing Tiantan Hospital, Capital Medical University, Beijing, China. ORCID: https://orcid.org/0000-0002-5960-4605. Tel: +86-10-59976611, Fax: +086-10-59978585, E-mail: abluelemon@163.com

Citation of this article: Xu N, Zhang L. Sepsis-associated Acute Kidney Injury: Epidemiology, Pathophysiology, and Management Strategies in Critical Care. *Nat Cell Sci* 2025;3(1):54–61. doi: 10.61474/ncs.2024.00011.

Abstract

Sepsis-associated acute kidney injury (SA-AKI) is a prevalent and serious complication of sepsis, accounting for over half of all sepsis-related complications and significantly contributing to morbidity and long-term mortality in intensive care units. This review elucidates the global epidemiology of SA-AKI, with an estimated incidence of approximately six million cases annually, and underscores diagnostic challenges stemming from temporal uncertainties and reliance on serum creatinine and urine output metrics. Our findings highlight that SA-AKI pathophysiology extends beyond traditional ischemia/reperfusion frameworks, incorporating inflammatory responses, microcirculatory dysfunction, and metabolic changes in renal tubules. Effective management requires a dual approach that addresses both sepsis control and targeted renal support. While continuous renal replacement therapy is preferred for hemodynamically unstable patients, the efficacy of various blood purification strategies requires further validation. Key prognostic indicators include injury severity, intervention timing, and renal recovery potential, with early detection and prompt intervention offering a means to improve outcomes. Emerging research on novel biomarkers holds promise for advancing diagnosis and prognosis, but optimal treatment strategies—particularly concerning the timing of renal replacement therapy initiation—require further refinement through future randomized controlled trials. This review highlights the urgent need for continued research to resolve these challenges and improve clinical outcomes for SA-AKI patients.

Keywords: Sepsis-associated acute kidney injury; SA-AKI; Epidemiology; Pathophysiology; Renal replacement therapy; RRT; Prognosis.

Introduction

Sepsis represents a clinical syndrome triggered by a dysregulated inflammatory response to infection, resulting in physiological and organ dysfunctions, 1,2 often leading to multi-organ failure. Acute kidney injury (AKI) is a critical health concern due to its high incidence, multifactorial etiology, and significant impact on morbidity and mortality, warranting focused research and clinical management efforts worldwide. AKI is the most common complication of sepsis, occurring in over 50% of cases.1 Among these, 15-20% of patients with Sepsis-Associated AKI (SA-AKI) require renal replacement therapy (RRT).^{2,3} Beyond its association with short-term mortality, SA-AKI also contributes to the development and progression of chronic kidney disease (CKD), end-stage renal disease, and increased long-term mortality.^{2,4} Management of SA-AKI primarily includes early fluid resuscitation, timely administration of appropriate antibiotics, vasopressors, infection source control, and supportive treatment for renal dysfunction.5 The Kidney Disease: Improv-

ing Global Outcomes (hereinafter referred to as KDIGO) AKI guidelines recommend continuous RRT (CRRT) for patients with hemodynamic instability.6 Moreover, large-scale intensive care unit (ICU) clinical investigations show that CRRT has become the preferred treatment for AKI patients with unstable hemodynamics, regardless of sepsis status.⁷ The pathophysiological mechanism of SA-AKI is incompletely understood. However, recent research suggests that dysregulated immune responses in sepsis trigger a cascade of inflammatory reactions, leading to a "cytokine storm" that accelerates SA-AKI progression.8 In addition to conventional CRRT, various novel extracorporeal blood purification therapies targeting endotoxins and cytokines have emerged. These approaches aim to regulate immune-inflammatory responses by selectively removing circulating endotoxins and cytokines. Current extracorporeal blood purification techniques for SA-AKI include CRRT, high cut-off filters, hemoperfusion/hemoadsorption, highly adsorptive hemofiltration, and continuous plasma filtration adsorption.9

Received: January 24, 2024 | Revised: May 08, 2024 | Accepted: March 21, 2025 | Published online: March 30, 2025



The Acute Dialysis Quality Initiative consensus distinguishes "sepsis-associated acute kidney injury" from "sepsis-induced acute kidney injury (SI-AKI)" to better characterize renal complications arising from sepsis. 10 SA-AKI broadly encompasses cases where AKI occurs in the setting of sepsis, reflecting the interplay between infection, systemic inflammation, and renal dysfunction. In contrast, SI-AKI specifically refers to kidney injury predominantly caused by sepsis-related pathophysiological derangements, such as systemic inflammatory responses and microcirculatory dysfunction. 11 Understanding these distinctions is crucial, given the heterogeneous nature of sepsis, which manifests in various phenotypes influenced by patient-specific factors such as genetic predisposition, comorbidities, and pathogen characteristics. 12 This diversity underscores the need for precision medicine approaches, which aim to tailor interventions based on individual sepsis phenotypes and their unique pathophysiological mechanisms, thereby optimizing treatment efficacy and improving patient outcomes.3 Precision medicine in SA-AKI and SI-AKI is expected to refine therapeutic strategies and improve management of these complex conditions by addressing the specific needs of diverse patient populations.

Epidemiology of SA-AKI

To date, very little is known about the epidemiology of SA-AKI. 13 Clinically, about 30% of sepsis patients develop AKI, with a global incidence of SA-AKI estimated at around six million cases per year, or nearly 0.1%. In the ICU, almost half of the patients with AKI develop sepsis. 9 A prospective cohort study involving 1,177 sepsis patients from multiple ICUs in Europe showed an AKI incidence of 51% and a mortality rate of 41%.¹⁴ A multicenter retrospective study found the incidence of AKI in sepsis patients to be 47.1%. This study included 1,243 patients with septic shock, of whom 50.4% were admitted to the emergency room due to AKI, and another 18.7% developed AKI within seven days. 15 Clinically, 70% of AKI patients are classified as Stage 2 or 3. AKI is also prevalent among patients without severe sepsis or septic shock. Approximately 34% of non-severe community-acquired pneumonia patients develop AKI. 16 These epidemiological data suggest a gradually increasing trend in the incidence of SA-AKI, which warrants greater research attention.17

Causes of SA-AKI

AKI in the context of SA-AKI represents a critical area of research due to its significant impact on patient morbidity and mortality rates in intensive care settings. SA-AKI is a common complication in septic patients, characterized by a sudden decline in renal function, which contributes to prolonged hospital stays, increased healthcare costs, and higher mortality rates. ¹⁸ Understanding the pathophysiological mechanisms and risk factors associated with SA-AKI is crucial for developing targeted interventions, improving patient outcomes, and enhancing critical care protocols. Current research underscores the complex interplay between systemic inflammation, hemodynamic instability, and cellular injury in the progression of SA-AKI, emphasizing the need

for early detection and personalized treatment strategies. ¹⁹ Addressing the causes and consequences of SA-AKI not only aids in the development of preventative measures but also enhances the overall management of sepsis, contributing to broader healthcare objectives of reducing the global burden of critical illnesses.

Diagnosis of SA-AKI

The exact onset of sepsis and AKI is often difficult to pinpoint. Both sepsis and its treatments can cause kidney damage.3 If AKI is directly caused by sepsis, it is termed sepsis-induced AKI, which is usually considered a subphenotype of SA-AKI. However, the interaction between sepsis and AKI is complex, with sepsis potentially causing AKI and critically ill patients possibly developing sepsis following AKI—suggesting that AKI may increase the risk of sepsis.²⁰ As such, the causality of SA-AKI remains to be fully determined. Moreover, diagnosing AKI based solely on serum creatinine and urine output has its limitations: sepsis can reduce peripheral blood perfusion, decreasing muscle blood supply, which reduces creatine production and delays an increase in serum creatinine levels, hindering early AKI detection. Inadequate urine output can also be diluted by aggressive fluid resuscitation, preventing timely diagnosis of AKI; the use of diuretics may further interfere with AKI diagnosis.²¹ Urine microscopy is a routine method for detecting kidney diseases. Compared to AKI patients caused by other factors, urine microscopy scores are higher in SA-AKI patients. When the urine microscopy score is ≥3, its specificity reaches up to 95%, but its sensitivity for detecting deteriorating AKI is relatively low (67%).6 This suggests that urinalysis may help determine the cause of AKI and provide prognostic information, but its sensitivity in diagnosing AKI and worsening AKI is not high.

Traditional biomarkers, such as serum creatinine, have limitations in the early diagnosis of acute kidney injury. In recent years, novel biomarkers like TIMP2IGFBP7 have been proposed for early AKI prediction and diagnosis. TIM-P2IGFBP7 assesses the risk of kidney injury by detecting the product of TIMP2 and IGFBP7 levels in urine, demonstrating strong predictive performance across various clinical settings. ^{22,23} As a cell cycle arrest biomarker, the TIMP2IGFBP7 concentration product in urine has been extensively studied for early AKI diagnosis. Studies indicate that TIMP2IGFBP7 has high sensitivity and specificity for predicting AKI in diverse clinical contexts. For instance, in a prospective study of 170 patients, the area under the curve for TIMP2IGFBP7 was 0.82, reflecting its robust predictive capability.²⁴ Another study found that TIMP2IGFBP7 exhibited high sensitivity and specificity for predicting AKI following cardiac surgery, with a threshold value of 0.265 (ng/ mL) □2/1,000.25 In addition to TIMP2IGFBP7, other novel biomarkers have also shown potential in early AKI diagnosis. Urinary biomarkers such as neutrophil gelatinase-associated lipocalin, kidney injury molecule-1, and liver-type fatty acid-binding protein have demonstrated good performance in early AKI prediction.²³ The combined use of these biomarkers can further enhance the accuracy of AKI predictions. Despite promising results from TIMP2IGFBP7 and other novel biomarkers, their clinical application faces several challenges. Firstly, the performance of these biomark-




Fig. 1. Epidemiology and Prognosis of SA-AKI. Sepsis-associated acute kidney injury (SA-AKI) presents a significant disease burden, with approximately 6 million cases reported annually. Nearly half of sepsis patients develop SA-AKI, and among those in the intensive care unit (ICU), the incidence is even higher, affecting 0.1% of ICU patients. Moreover, 70% of ICU patients with SA-AKI meet stage 2–3 severity criteria according to the Kidney Disease: Improving Global Outcomes (KDIGO) guidelines. This flowchart illustrates the potential progression and outcomes of SA-AKI. Upon hospital admission, patients present with SA-AKI without a clearly defined sequence of preceding sepsis or AKI events. Following admission, SA-AKI patients may follow one of three pathways over time: (1) early recovery of kidney function, (2) late recovery of kidney function, or (3) no recovery of kidney function. The third pathway, where kidney function does not recover, may further progress to acute kidney disease (AKD) or chronic kidney disease (CKD) or result in death. The flowchart depicts these three potential trajectories for SA-AKI patients, starting from the point of hospital admission when AKI is first identified.

ers can vary depending on the clinical setting and patient population. For instance, while TIMP2*IGFBP7 performed well in predicting AKI post-cardiac surgery in some studies, its effectiveness may vary in others. ²⁶ Secondly, detecting these biomarkers requires specific equipment and techniques, limiting their application in resource-constrained settings. Furthermore, the critical thresholds and optimal timing for biomarker testing require additional research and validation. ²³

Course of SA-AKI

The precise onset of SA-AKI is currently unclear. In patients with sepsis, the existence of AKI should be suspected; conversely, the presence of sepsis should be considered in AKI patients (Fig. 1). AKI can occur simultaneously with sepsis at admission or manifest during hospitalization. Even when AKI is not present at the time of admission in septic patients, optimal resuscitation and appropriate sepsis treatment can still act preventatively against AKI.21 Once diagnosed with SA-AKI, there should be an emphasis on thorough monitoring and timely organ-supportive treatment to prevent further renal injury. SA-AKI may reverse within the first week of admission, indicating a generally positive prognosis.27 However, some patients may experience one or more episodes of AKI, highlighting the importance of continuous monitoring and avoidance of nephrotoxic injury, even if initial AKI has been reversed or recovered. Even patients who fully recover from SA-AKI remain at risk of developing CKD and other outcomes, including recurrent sepsis. 28 Patients not fully recovered from AKI within seven days will be categorized as having acute kidney disease, which may later recover or progress to CKD and is associated with a poor long-term prognosis.²⁹ Post-discharge survivors of SA-AKI should be followed up by nephrologists for continued management to monitor the occurrence of CKD and other adverse outcomes.

Pathophysiological mechanism of SA-AKI

The pathophysiological mechanism of SA-AKI is still not completely understood. It is generally believed that SA-AKI is due to a reduction in bilateral renal glomerular blood perfusion, resulting in acute renal tubular ischemia and even apoptosis of renal tubular epithelial cells (RTEC). Additionally, inflammatory responses, metabolic adaptations, and microcirculatory dysfunction play significant roles in the organ damage phase of sepsis.²⁹

Ischemia/reperfusion (I/R) injury

Hypoperfusion and shock are primary causes of AKI, with I/R injury leading to acute tubular necrosis and extensive cell death. However, during SA-AKI, various mechanisms operate beyond I/R injury.³⁰ For instance, SA-AKI may occur even when renal hemodynamics are stable.

Inflammatory response

The inflammatory response is a defense mechanism that forms after a pathogen invades the host. During sepsis, inflammatory mediators such as pathogen-associated molecular patterns and damage-associated molecular patterns are released into the bloodstream.³¹ These molecules bind to pattern recognition receptors, such as Toll-like receptors (TLRs), on immune cells, triggering a downstream signal-

ing cascade that leads to the synthesis and release of proinflammatory molecules. RTECs also express TLRs, specifically TLR2 and TLR4. 32,33 Upon pathogen-associated molecular pattern activation, proximal RTECs exhibit increased oxidative stress, produce reactive oxygen species, and cause mitochondrial damage. There is evidence that RTECs initiate signal transduction through paracrine secretion. 34 Furthermore, pathological observations revealed increased monocyte infiltration in the glomeruli and peritoneal regions of sepsis animal models, confirming the role of inflammation in AKI. 35

Microcirculatory dysfunction

Studies have shown that even with stable renal blood supply, sepsis patients can still experience alterations in renal microcirculation.³⁶ Features of SA-AKI include changes in the homogeneity of microcirculatory blood flow, manifested as a reduced proportion of capillary blood flow and an increased proportion of stalled flow.³⁷ Multiple mechanisms may cause microcirculatory dysfunction, such as sympathetic and parasympathetic nervous system responses, endothelial injury, sloughing of the glycocalyx, and activation of the coagulation cascade. Endothelial injury reinforces leukocyte infiltration and platelet adhesion, while slowing blood flow, leading to capillary occlusion, microthrombus formation, and prolonged exposure of RTEC to activated inflammatory mediators. The increase in vascular permeability and impaired dilation function related to endothelial injury results in perivascular interstitial edema.³⁸ This alters perfusion to the RTEC by increasing the oxygen diffusion distance from the capillaries to the RTEC and elevating the output pressure of the vein. Changes in microcirculatory hemodynamics may also play a key role in SA-AKI. The glomerular filtration rate is determined by intraglomerular hydrostatic pressure and is independent of renal blood flow changes; hence, the mechanism of afferent arteriolar constriction and efferent arteriolar dilation is suggested to explain the decrease in glomerular pressure, which results in reduced glomerular filtration rate.²⁸

Additionally, during AKI, renal blood flow undergoes redistribution, diverting away from the medulla, bypassing the glomerulus, and directly connecting the incoming substances to the efferent capillaries of the arteriole. This partially explains blood shunting during SA-AKI. However, it remains unclear how or when these shunting pathways open. In summary, the increase in shunting and redistribution of blood flow explains the possible presence of mechanisms beyond ischemia during SA-AKI.³⁹

Metabolic adaptation

During SA-AKI, metabolic adaptation can occur, where a method is adopted that maintains the minimum function of cells and organs to ensure cell survival. Various theories have been proposed to explain the metabolic adaptation of RTEC cells during sepsis, most of which are mitochondria-mediated and characterized by optimized energy consumption. Inflammation is associated with optimized energy consumption, implying a reduction in the efficiency of energy utilization in protein synthesis or ion transport, thereby maintaining only energy use for crucial cell functions and simultaneously avoiding cell death.⁴⁰ During inflammation, the expression of renal tubule transporter proteins is down-

regulated, and renal tubule ion transport is reduced, suggesting that metabolic adaptation can re-determine priority in energy consumption, which is an adaptive mechanism. ⁴¹ The exact mechanisms of how metabolic adaptation occurs are still unclear.

Lipid metabolism alterations

The alterations in lipid metabolism observed in sepsis are complex and multifaceted, affecting not only the progression of the disease but also potentially providing new therapeutic targets. 42 Sepsis is a systemic inflammatory response syndrome triggered by infection, characterized by an excessive host response that leads to multi-organ dysfunction. In patients with sepsis, there are significant changes in lipid metabolism, particularly a decrease in high-density lipoprotein (HDL) levels, which are closely linked to the exacerbation of the inflammatory response and endothelial dysfunction.⁴³ The role of HDL in sepsis extends beyond cholesterol transport, encompassing anti-inflammatory, anti-apoptotic, antithrombotic, and anti-infective properties. During sepsis, both the serum levels and composition of HDL are significantly altered, potentially leading to dysregulation of host immune responses and endothelial dysfunction. Research indicates that the anti-inflammatory capacity of HDL is suppressed in sepsis, possibly due to the release of inflammatory mediators and increased oxidative stress. Furthermore, changes in lipid metabolism may influence the prognosis of sepsis. For instance, low HDL levels are associated with poor outcomes in sepsis patients, suggesting HDL's potential as a biomarker for predicting sepsis severity and prognosis. Additionally, some studies have explored the potential of treating sepsis by modulating lipid metabolism.42,43

Treatment of SA-AKI

Control of the pathogenesis and early, appropriate use of antimicrobial drugs remain vital in the treatment of sepsis and can also prevent further renal injury. Delay in the administration of antimicrobial drugs is associated with the early development of AKI. However, caution should be exercised when using certain nephrotoxic drugs involved in treatment, such as aminoglycosides and vancomycin. The KDIGO guideline suggests that diagnostic drugs such as amphotericin B and intravenous radiographic contrast agents should be used cautiously to avoid renal injury, with strict drug monitoring considered when necessary. A bundled treatment approach for SA-AKI may be beneficial, including infection source control, appropriate antibiotics, fluid therapy, vasopressors, avoidance of nephrotoxic insults, and RRT (Fig. 2).

Fluid resuscitation

The use of vasoactive medications after fluid resuscitation forms the foundation of shock treatment. Three landmark clinical trials conducted among septic shock patients have consistently indicated that resuscitation strategies do not confer advantages in terms of reducing mortality rates or the need for RRT.⁵ Additionally, the influence of early goal-directed therapy, alternative resuscitation strategies, or conventional therapy on the development of new AKI, the severity of AKI, fluid overload, the need for RRT, or renal recovery is minimal.⁴⁴




Fig. 2. Treatment of SA-AKI. This diagram outlines a bundled management approach for sepsis-associated acute kidney injury (SA-AKI). The bundle consists of six key interventions: (1) infection source control, (2) appropriate antibiotic use, (3) fluid therapy, (4) vasopressor administration, (5) avoidance of nephrotoxic insults, and (6) renal replacement therapy (RRT).

Types of resuscitation fluids

Every type of intravenous fluid infusion has the potential to lead to fluid overload and renal burden. A prospective observational study indicated that patients with AKI had a greater cumulative fluid intake in the first three days, suggesting that fluid overload can function as an independent risk factor for the severity of AKI. Evidence demonstrates that the use of hydroxyethyl starch elevates the incidence of AKI, the need for RRT, and the risk of mortality, especially when administered in large quantities to patients with sepsis. 14 In contrast to crystalloids, gelatin may correlate with a higher occurrence rate of AKI and a greater requirement for RRT among septic patients. 45 Albumin solution is usually considered safe for sepsis resuscitation; however, recent observational research involving shock patients revealed a significant correlation between early, substantial albumin infusion postsurgery and the development of AKI.46

An observational study involving ICU patients, including those experiencing septic shock, demonstrated a reduced incidence and mortality rate of AKI when using Ringer's solution. Notably, even within sepsis subgroups, there was no statistically significant difference in AKI incidence or the proportion requiring RRT when comparing patients treated with normal saline to those treated with balanced crystalloid solutions.⁴⁷ Other reports indicate that the risk of major adverse renal events is reduced among patients treated with balanced crystalloid solutions.⁴⁸ Therefore, balanced solutions should be used instead of normal saline, particularly in patients with sepsis.⁴⁹

Vasoactive drugs

Norepinephrine is the drug of choice for the treatment of septic shock. Dopamine, compared to norepinephrine, is associated with more adverse events and is not recommended for renal protection. 50 Vasopressin is associated with a lower rate of RRT and does not appear to increase the risk of AKI. 51 A large randomized controlled trial for septic shock showed that a blood pressure level of 80–85 mmHg (1 mmHg \approx 0.133 kPa) reduced the need for RRT in patients compared to a target mean arterial pressure of 65–70 mmHg, but no survival benefit was observed. 52

Initiation timing of renal replacement therapy

Early application of RRT has been shown to significantly benefit patients with severe AKI within 90 days of initiating RRT.46 Concurrently, early treatment substantially reduces the incidence of major adverse kidney events and mortality within a year, contributing to the restoration of kidney function. 53 However, the findings of a trial comparing the effectiveness of early versus delayed dialysis in patients with septic shock and severe AKI in the ICU demonstrated no statistically significant difference in mortality rates between these two RRT approaches.⁵⁴ Alarmingly, 9% of patients died regardless of whether early or delayed initiation was used.53 The Standard versus Accelerated Initiation of Renal-Replacement Therapy in Acute Kidney Injury trial investigated the optimal timing for initiating kidney replacement therapy in patients with acute kidney injury. The trial compared two strategies: accelerated initiation versus standard initiation of KRT. The findings revealed a neutral outcome concerning 90-day all-cause mortality, indicating that early initiation of KRT does not confer additional survival benefits compared to the standard initiation strategy in critically ill patients.⁵⁵ The appropriate timing for the initiation of RRT is still under debate.

Blood purification

Clinical evidence on blood purification remains limited, with most studies not having measured target solutes. It is still unclear whether inflammatory mediators can be removed, and even if they can, whether their removal benefits the prognosis or survival. In patients with septic shock, cytokine levels are highly variable, and high endogenous clearance rates make their levels dynamic. 56 The EUPHRATES trial evaluated the efficacy of using a blood purification device in patients with severe sepsis and septic shock. While the primary outcome did not demonstrate a significant improvement in overall survival rates, potential benefits were observed in certain subgroups of patients.⁵⁷ Furthermore, the development of novel therapies based on adsorption devices has drawn attention to the clearance of specific mediators involved in organ dysfunction and antibiotic removal. The joint commission of the Italian Society of Anesthesiology and Critical Care and the Italian Society of Nephrology has addressed some of these issues, proposed clinical practice recommendations, and developed a framework for future research in this field.58

High-volume hemofiltration (HVHF), defined as a convective dose of >35 mL·kg⁻¹·h⁻¹, has been explored in several studies regarding its efficacy in treating SA-AKI patients, with varying results in mortality rates. A multi-center randomized controlled trial involving 140 SA-AKI patients, divided into two groups receiving either HVHF at a rate of 70 mL·kg⁻¹·h⁻¹ or standard hemofiltration at 35 mL·kg⁻¹·h⁻¹ for 96 h, demonstrated no significant statistical difference in patient mortality at 28, 60, and 90 days between the two treatment methods.⁵⁹ Another randomized controlled trial comparing the treatment effects of varying dosages of HVHF in SA-AKI patients also showed no significant statistical difference in mortality or renal outcomes between patients treated with 85 mL·kg⁻¹·h⁻¹ and 50 mL·kg⁻¹·h⁻¹ HVHF.⁶⁰ Additionally, a meta-analysis incorporating four HVHF trials involving SA-AKI patients indicated no benefit in 28-day survival rates and an increased incidence of hypophosphatemia and hypokalemia.61

Polymyxin B hemoperfusion (PMX-HP), utilized as an adjunctive measure in sepsis treatment, can remove endotoxins from the blood circulation. Various randomized controlled trials have confirmed the efficacy of PMX-HP in sepsis, though some studies have yielded contrary findings. A European multi-center pilot study involving 36 intra-abdominal sepsis surgical patients demonstrated that a 2-h PMX-HP treatment could only improve left ventricular function and decrease the need for RRT.62,63 Research that included 64 severe intra-abdominal sepsis patients showed benefits of PMX-HP treatment in terms of hemodynamics, organ function, and 28-day survival rate.64 However, a subsequent larger randomized controlled trial yielded negative results.65 Additionally, the immunomodulatory role of PMX-HP was reported in another randomized controlled trial, validating its ability to enhance the recruitment of monocytes and neutrophils in sepsis patients, but without any benefit in reducing

mortality or improving renal outcomes.66

Prognosis of SA-AKI

SA-AKI is closely associated with poor clinical outcomes, including acute kidney disease. CKD, or death (Fig. 1). In critically ill patients with AKI, those with SA-AKI have a higher risk of in-hospital death and a longer hospital stay. 9,67 In patients with SA-AKI that reverses within 24 h post-shock, there is a reduced in-hospital mortality rate. The long-term prognosis of SA-AKI patients depends on the primary disease, the severity of AKI, and the baseline condition at discharge. Patients who recover from AKI have significantly improved survival rates, but there remains a possibility of CKD. end-stage renal disease, and death. In a follow-up of over a year for 105 survivors, it was found that the rates of CKD in patients with reversed AKI, recovered AKI, and unrecovered AKI were 21%, 30%, and 79%, respectively.68 Another longterm follow-up study of 19,276 survivors showed that the 15year mortality rates of suspected SA-AKI by renal recovery status at discharge were 57.1% for no recovery, 47.9% for recovery, and 36.5% for no AKI.54 The severity of AKI, timely treatment, and recovery status during hospitalization are critical factors in determining the prognosis.69

Conclusions

In summary, one-third of sepsis patients develop AKI, with half of these cases being severe. Close monitoring and timely treatment are crucial for SA-AKI to prevent further kidney injury and improve prognosis. The main mechanisms include reduced renal perfusion, inflammation, and impaired microcirculation. Effective sepsis treatment involves addressing the cause and using antibiotics, fluid resuscitation, and vasopressors appropriately. Research is advancing in finding biomarkers for early SA-AKI diagnosis, such as neutrophil gelatinase-associated lipocalin and kidney injury molecule-1, which may improve future diagnostics. Further studies, particularly on the timing of RRT, are essential to enhance treatment outcomes.

Acknowledgments

We acknowledge Figdraw (www.figdraw.com) for the material provided for the picture production.

Funding

None.

Conflict of interest

One of the authors, Linlin Zhang, serves as an editorial board member of *Nature Cell and Science*. The authors declare no other conflicts of interest.

Author contributions

Conceptualization: L.Z.; Project administration: L.Z; Resources: N.X.; Visualization: N.X.; Writing – original draft: N.X.;

Writing - review & editing: L.Z.

References

- Kellum JA, Formeck CL, Kernan KF, Gómez H, Carcillo JA. Subtypes and Mimics of Sepsis. Crit Care Clin 2022;38(2):195–211. doi:10.1016/j. ccc.2021.11.013, PMID:35369943.
- [2] Manrique-Caballero CL, Del Rio-Pertuz G, Gomez H. Sepsis-Associated Acute Kidney Injury. Crit Care Clin 2021;37(2):279–301. doi:10.1016/j. ccc.2020.11.010, PMID:33752856.
- [3] Zarbock A, Nadim MK, Pickkers P, Gomez H, Bell S, Joannidis M, et al. Sepsis-associated acute kidney injury: consensus report of the 28th Acute Disease Quality Initiative workgroup. Nat Rev Nephrol 2023;19(6):401–417. doi:10.1038/s41581-023-00683-3, PMID:368 23168.
- [4] Yue S, Li S, Huang X, Liu J, Hou X, Zhao Y, et al. Machine learning for the prediction of acute kidney injury in patients with sepsis. J Transl Med 2022;20(1):215. doi:10.1186/s12967-022-03364-0, PMID:35562803.
- [5] Peerapornratana S, Manrique-Caballero CL, Gómez H, Kellum JA. Acute kidney injury from sepsis: current concepts, epidemiology, pathophysiology, prevention and treatment. Kidney Int 2019;96(5):1083–1099. doi:10.1016/j.kint.2019.05.026, PMID:31443997.
- [6] White KC, Serpa-Neto A, Hurford R, Clement P, Laupland KB, See E, et al. Sepsis-associated acute kidney injury in the intensive care unit: incidence, patient characteristics, timing, trajectory, treatment, and associated outcomes. A multicenter, observational study. Intensive Care Med 2023;49(9):1079–1089. doi:10.1007/s00134-023-07138-0, PMID:37432520.
- [7] Yoshimoto K, Komaru Y, Iwagami M, Doi K. Acute Kidney Injury in Sepsis: Evidence From Asia. Semin Nephrol 2020;40(5):489–497. doi:10.1016/j.semnephrol.2020.08.005, PMID:33334462.
- [8] Wang K, Xie S, Xiao K, Yan P, He W, Xie L. Biomarkers of Sepsis-Induced Acute Kidney Injury. Biomed Res Int 2018;2018:6937947. doi:10.1155/2018/6937947, PMID:29854781.
- [9] Kellum JA, Wen X, de Caestecker MP, Hukriede NA. Sepsis-Associated Acute Kidney Injury: A Problem Deserving of New Solutions. Nephron 2019;143(3):174–178. doi:10.1159/000500167, PMID:31018211.
- [10] Sun S, Chen R, Dou X, Dai M, Long J, Wu Y, et al. Immunoregulatory mechanism of acute kidney injury in sepsis: A Narrative Review. Biomed Pharmacother 2023;159:114202. doi:10.1016/j.biopha.2022.114202, PMID:36621143.
- [11] Sahin Aktura S, Sahin K, Tumkaya L, Mercantepe T, Topcu A, Pinarbas E, et al. The Nephroprotective Effect of Punica granatum Peel Extract on LPS-Induced Acute Kidney Injury. Life (Basel) 2024;14(10):1316. doi:10.3390/life14101316, PMID:39459616.
- [12] Kounatidis D, Tzivaki I, Daskalopoulou S, Daskou A, Adamou A, Rigatou A, et al. Sepsis-Associated Acute Kidney Injury: What's New Regarding Its Diagnostics and Therapeutics? Diagnostics (Basel) 2024;14(24):2845. doi:10.3390/diagnostics14242845, PMID:39767206.
- [13] Wang H, Ji X, Wang AY, Wu PK, Liu Z, Dong L, et al. Epidemiology of Sepsis-Associated Acute Kidney Injury in Beijing, China: A Descriptive Analysis. Int J Gen Med 2021;14:5631–5649. doi:10.2147/IJGM. S320768, PMID:34548815.
- [14] Chen WY, Cai LH, Zhang ZH, Tao LL, Wen YC, Li ZB, et al. The timing of continuous renal replacement therapy initiation in sepsis-associated acute kidney injury in the intensive care unit: the CRTSAKI Study (Continuous RRT Timing in Sepsis-associated AKI in ICU): study protocol for a multicentre, randomised controlled trial. BMJ Open 2021;11(2):e040718. doi:10.1136/bmjopen-2020-040718, PMID:33608398.
- [15] Xu X, Nie S, Liu Z, Chen C, Xu G, Zha Y, et al. Epidemiology and Clinical Correlates of AKI in Chinese Hospitalized Adults. Clin J Am Soc Nephrol 2015;10(9):1510–1518. doi:10.2215/CJN.02140215, PMID:26231194.
- [16] Hoste EAJ, Kellum JA, Selby NM, Zarbock A, Palevsky PM, Bagshaw SM, et al. Global epidemiology and outcomes of acute kidney injury. Nat Rev Nephrol 2018;14(10):607–625. doi:10.1038/s41581-018-0052-0, PMID:30135570.
- [17] Rossaint J, Zarbock A. Acute kidney injury: definition, diagnosis and epidemiology. Minerva Urol Nefrol 2016;68(1):49–57. PMID:26364570.

- [18] Yang L, Xing G, Wang L, Wu Y, Li S, Xu G, et al. Acute kidney injury in China: a cross-sectional survey. Lancet 2015;386(10002):1465–1471. doi:10.1016/S0140-6736(15)00344-X, PMID:26466051.
- [19] Aguilar MG, AlHussen HA, Gandhi PD, Kaur P, Pothacamuri MA, Talikoti MAH, et al. Sepsis-Associated Acute Kidney Injury: Pathophysiology and Treatment Modalities. Cureus 2024;16(12):e75992. doi:10.7759/ cureus.75992, PMID:39834999.
- [20] Keir I, Kellum JA. Acute kidney injury in severe sepsis: pathophysiology, diagnosis, and treatment recommendations. J Vet Emerg Crit Care (San Antonio) 2015;25(2):200–209. doi:10.1111/vec.12297, PMID:25845505.
- [21] Molinari L, Del Rio-Pertuz G, Smith A, Landsittel DP, Singbartl K, Palevsky PM, et al. Utility of Biomarkers for Sepsis-Associated Acute Kidney Injury Staging. JAMA Netw Open 2022;5(5):e2212709. doi:10.1001/ jamanetworkopen.2022.12709, PMID:35583867.
- [22] Fiorentino M, Castellano G, Kellum JA. Differences in acute kidney injury ascertainment for clinical and preclinical studies. Nephrol Dial Transplant 2017;32(11):1789–1805. doi:10.1093/ndt/gfx002, PMID:28371878.
- [23] Cheng L, Jia HM, Zheng X, Jiang YJ, Zhang TE, Li WX. Urinary cell cycle biomarkers for the prediction of renal non-recovery in patients with septic acute kidney injury: a prospective study. Clin Exp Nephrol 2023;27(12):1051–1059. doi:10.1007/s10157-023-02397-z, PMID: 37656396.
- [24] Inotani S, Kashio T, Osakabe Y, Matsumoto T, Nagao Y, Ishihara M, et al. Efficacy of urinary [TIMP-2]@[IGFBP7], L-FABP, and NGAL levels for predicting community-acquired acute kidney injury in Japanese patients: a single-center, prospective cohort study. Clin Exp Nephrol 2025. doi:10.1007/s10157-025-02641-8, PMID:39982587.
- [25] Xie Y, Guo Q, Yang B, Xie P, Zhang J, Lu W, et al. Tissue Inhibitor Metalloproteinase-2-IGF-Binding Protein 7 for the Prediction of Acute Kidney Injury following Cardiac Surgery. Cardiorenal Med 2024;14(1):251– 260. doi:10.1159/000538031, PMID:38588650.
- [26] Yu JT, Hu XW, Yang Q, Shan RR, Zhang Y, Dong ZH, et al. Insulin-like growth factor binding protein 7 promotes acute kidney injury by alleviating poly ADP ribose polymerase 1 degradation. Kidney Int 2022;102(4):828–844. doi:10.1016/j.kint.2022.05.026, PMID:35752325.
- [27] Nusshag C, Wei C, Hahm E, Hayek SS, Li J, Samelko B, *et al.* suPAR links a dysregulated immune response to tissue inflammation and sepsis-induced acute kidney injury. JCI Insight 2023;8(7):e165740. doi:10.1172/jci.insight.165740, PMID:37036003.
- [28] Bouglé A, Duranteau J. Pathophysiology of sepsis-induced acute kidney injury: the role of global renal blood flow and renal vascular resistance. Contrib Nephrol 2011;174:89–97. doi:10.1159/000329243, PMID:21921613.
- [29] Ozkaya PY, Taner S, Ersayoğlu I, Turan B, Yildirim Arslan S, Karapinar B, et al. Sepsis associated acute kidney injury in pediatric intensive care unit. Ther Apher Dial 2023;27(1):73–82. doi:10.1111/1744-9987.13928, PMID:36101487.
- [30] Hasegawa S, Inoue T, Nakamura Y, Fukaya D, Uni R, Wu CH, et al. Activation of Sympathetic Signaling in Macrophages Blocks Systemic Inflammation and Protects against Renal Ischemia-Reperfusion Injury. J Am Soc Nephrol 2021;32(7):1599–1615. doi:10.1681/ASN.2020121723, PMID:33875568.
- [31] Zhu Y, Xu D, Deng F, Yan Y, Li J, Zhang C, et al. Angiotensin (1-7) Attenuates Sepsis-Induced Acute Kidney Injury by Regulating the NF-κB Pathway. Front Pharmacol 2021;12:601909. doi:10.3389/fphar.2021.601909, PMID:33746749.
- [32] Zhong Y, Wu S, Yang Y, Li GQ, Meng L, Zheng QY, et al. LIGHT aggravates sepsis-associated acute kidney injury via TLR4-MyD88-NF-κB pathway. J Cell Mol Med 2020;24(20):11936–11948. doi:10.1111/jcmm.15815, PMID:32881263.
- [33] Wang M, Wei J, Shang F, Zang K, Zhang P. Down-regulation of lncRNA SNHG5 relieves sepsis-induced acute kidney injury by regulating the miR-374a-3p/TLR4/NF-κB pathway. J Biochem 2021;169(5):575–583. doi:10.1093/jb/mvab008, PMID:33479745.
- [34] Fani F, Regolisti G, Delsante M, Cantaluppi V, Castellano G, Gesualdo L, et al. Recent advances in the pathogenetic mechanisms of sepsis-associated acute kidney injury. J Nephrol 2018;31(3):351–359. doi:10.1007/s40620-017-0452-4, PMID:29273917.

- [35] Keshari RS, Popescu NI, Silasi R, Regmi G, Lupu C, Simmons JH, et al. Complement C5 inhibition protects against hemolytic anemia and acute kidney injury in anthrax peptidoglycan-induced sepsis in baboons. Proc Natl Acad Sci U S A 2021;118(37):e2104347118. doi:10.1073/pnas.2104347118, PMID:34507997.
- [36] Kuwabara S, Goggins E, Okusa MD. The Pathophysiology of Sepsis-Associated AKI. Clin J Am Soc Nephrol 2022;17(7):1050–1069. doi:10.2215/CJN.00850122, PMID:35764395.
- [37] Ma S, Evans RG, Iguchi N, Tare M, Parkington HC, Bellomo R, et al. Sepsis-induced acute kidney injury: A disease of the microcirculation. Microcirculation 2019;26(2):e12483. doi:10.1111/micc.12483, PMID:29908046.
- [38] Prowle JR, Bellomo R. Sepsis-associated acute kidney injury: macrohemodynamic and microhemodynamic alterations in the renal circulation. Semin Nephrol 2015;35(1):64–74. doi:10.1016/j.semnephrol.2015.01.007, PMID:25795500.
- [39] Ergin B, Kapucu A, Demirci-Tansel C, Ince C. The renal microcirculation in sepsis. Nephrol Dial Transplant 2015;30(2):169–177. doi:10.1093/ ndt/gfu105, PMID:24848133.
- [40] Umbro I, Gentile G, Tinti F, Muiesan P, Mitterhofer AP. Recent advances in pathophysiology and biomarkers of sepsis-induced acute kidney injury. J Infect 2016;72(2):131–142. doi:10.1016/j.jinf.2015.11.008, PMID:26702738.
- [41] Rosenberger C, Rosen S, Heyman SN. Renal parenchymal oxygenation and hypoxia adaptation in acute kidney injury. Clin Exp Pharmacol Physiol 2006;33(10):980–988. doi:10.1111/j.1440-1681.2006. 04472.x, PMID:17002677.
- [42] Stasi A, Franzin R, Fiorentino M, Squiccimarro E, Castellano G, Gesualdo L. Multifaced Roles of HDL in Sepsis and SARS-CoV-2 Infection: Renal Implications. Int J Mol Sci 2021;22(11):5980. doi:10.3390/ijms 22115980, PMID:34205975.
- [43] Stasi A, Fiorentino M, Franzin R, Staffieri F, Carparelli S, Losapio R, et al. Beneficial effects of recombinant CER-001 high-density lipoprotein infusion in sepsis: results from a bench to bedside translational research project. BMC Med 2023;21(1):392. doi:10.1186/s12916-023-03057-5, PMID:37915050.
- [44] Zarbock A, Koyner JL, Gomez H, Pickkers P, Forni L, Acute Disease Quality Initiative group. Sepsis-associated acute kidney injury-treatment standard. Nephrol Dial Transplant 2023;39(1):26–35. doi:10.1093/ndt/gfad142, PMID:37401137.
- [45] Karkar A. Continuous renal replacement therapy: Principles, modalities, and prescription. Saudi J Kidney Dis Transpl 2019;30(6):1201–1209. doi:10.4103/1319-2442.275463, PMID:31929266.
- [46] Li X, Liu C, Mao Z, Li Q, Zhou F. Timing of renal replacement therapy initiation for acute kidney injury in critically ill patients: a systematic review of randomized clinical trials with meta-analysis and trial sequential analysis. Crit Care 2021;25(1):15. doi:10.1186/s13054-020-03451-y, PMID:33407756.
- [47] Wu S, Xu T, Wu C, Lei X, Tian X. Continuous renal replacement therapy in sepsis-associated acute kidney injury: Effects on inflammatory mediators and coagulation function. Asian J Surg 2021;44(10):1254–1259. doi:10.1016/j.asjsur.2021.02.004, PMID:34088582.
- [48] Koyner JL. Sepsis and Kidney Injury. Contrib Nephrol 2021;199:56–70. doi:10.1159/000517701, PMID:34343994.
- [49] Serpa Neto A, Veelo DP, Peireira VG, de Assunção MS, Manetta JA, Espósito DC, et al. Fluid resuscitation with hydroxyethyl starches in patients with sepsis is associated with an increased incidence of acute kidney injury and use of renal replacement therapy: a systematic review and meta-analysis of the literature. J Crit Care 2014;29(1):185. e1–185.e7. doi:10.1016/j.jcrc.2013.09.031, PMID:24262273.
- [50] Al-Husinat L, Alsabbah A, Hmaid AA, Athamneh R, Adwan M, Hourani MN, et al. Norepinephrine May Exacerbate Septic Acute Kidney Injury: A Narrative Review. J Clin Med 2023;12(4):1373. doi:10.3390/icm12041373. PMID:36835909.
- [51] Jozwiak M. Alternatives to norepinephrine in septic shock: Which agents and when? J Intensive Med 2022;2(4):223–232. doi:10.1016/j. jointm.2022.05.001, PMID:36788938.
- [52] Wei S, Gao Y, Dai X, Fu W, Cai S, Fang H, et al. SIRT1-mediated HMGB1 deacetylation suppresses sepsis-associated acute kidney injury. Am J Physiol Renal Physiol 2019;316(1):F20–F31. doi:10.1152/ajpre-nal.00119.2018, PMID:30379096.

- [53] Wang Z, Zhang L, Xu F, Han D, Lyu J. The association between continuous renal replacement therapy as treatment for sepsis-associated acute kidney injury and trend of lactate trajectory as risk factor of 28-day mortality in intensive care units. BMC Emerg Med 2022;22(1):32. doi:10.1186/s12873-022-00589-6, PMID:35227200.
- [54] Kim IY, Kim S, Ye BM, Kim MJ, Kim SR, Lee DW, et al. Procalcitonin decrease predicts survival and recovery from dialysis at 28 days in patients with sepsis-induced acute kidney injury receiving continuous renal replacement therapy. PLoS One 2022;17(12):e0279561. doi:10.1371/journal.pone.0279561, PMID:36574383.
- [55] Zampieri FG, da Costa BR, Vaara ST, Lamontagne F, Rochwerg B, Nichol AD, et al. A Bayesian reanalysis of the Standard versus Accelerated Initiation of Renal-Replacement Therapy in Acute Kidney Injury (STARRT-AKI) trial. Crit Care 2022;26(1):255. doi:10.1186/s13054-022-04120-y, PMID:36008827.
- [56] Monard C, Rimmelé T, Ronco C. Extracorporeal Blood Purification Therapies for Sepsis. Blood Purif 2019;47(Suppl 3):1–14. doi:10.1159/000499520, PMID:30974444.
- [57] Rachoin JS, Foster D, Dellinger RP. Endotoxin removal: how far from the evidence? From EUPHAS to EUPHRATES. Contrib Nephrol 2010;167:111–118. doi:10.1159/000315925, PMID:20519905.
- [58] De Rosa S, Marengo M, Fiorentino M, Fanelli V, Brienza N, Fiaccadori E, et al. Extracorporeal blood purification therapies for sepsis-associated acute kidney injury in critically ill patients: expert opinion from the SIAARTI-SIN joint commission. J Nephrol 2023;36(7):1731–1742. doi:10.1007/s40620-023-01637-5, PMID:37439963.
- [59] Bellomo R, Kellum JA, Ronco C, Wald R, Martensson J, Maiden M, et al. Acute kidney injury in sepsis. Intensive Care Med 2017;43(6):816–828. doi:10.1007/s00134-017-4755-7, PMID:28364303.
- [60] Meng SQ, Yang WB, Liu JG, Yuan JY, Zhang K, Ding WY, et al. Evaluation of the application of high volume hemofiltration in sepsis combined with acute kidney injury. Eur Rev Med Pharmacol Sci 2018;22(3):715– 720. doi:10.26355/eurrev_201802_14298, PMID:29461600.
- [61] Clark E, Molnar AO, Joannes-Boyau O, Honoré PM, Sikora L, Bagshaw SM. High-volume hemofiltration for septic acute kidney injury: a systematic review and meta-analysis. Crit Care 2014;18(1):R7. doi:10.1186/cc13184, PMID:24398168.
- [62] Lankadeva YR, Okazaki N, Evans RG, Bellomo R, May CN. Renal Medullary Hypoxia: A New Therapeutic Target for Septic Acute Kidney Injury? Semin Nephrol 2019;39(6):543–553. doi:10.1016/j.semnephrol.2019.10.004. PMID:31836037.
- [63] Cantaluppi V, Assenzio B, Pasero D, Romanazzi GM, Pacitti A, Lanfranco G, et al. Polymyxin-B hemoperfusion inactivates circulating proapoptotic factors. Intensive Care Med 2008;34(9):1638–1645. doi:10.1007/s00134-008-1124-6, PMID:18463848.
- [64] Wan L, Bagshaw SM, Langenberg C, Saotome T, May C, Bellomo R. Pathophysiology of septic acute kidney injury: what do we really know? Crit Care Med 2008;36(4 Suppl):S198–S203. doi:10.1097/ CCM.0b013e318168ccd5, PMID:18382194.
- [65] Jacobs R, Honore PM, Joannes-Boyau O, Boer W, De Regt J, De Waele E, et al. Septic acute kidney injury: the culprit is inflammatory apoptosis rather than ischemic necrosis. Blood Purif 2011;32(4):262–265. doi:10.1159/000330244, PMID:21860231.
- [66] Chen JJ, Lai PC, Lee TH, Huang YT. Blood Purification for Adult Patients With Severe Infection or Sepsis/Septic Shock: A Network Meta-Analysis of Randomized Controlled Trials. Crit Care Med 2023;51(12):1777– 1789. doi:10.1097/CCM.00000000000005991. PMID:37470680.
- [67] Liang M, Ren X, Huang D, Ruan Z, Chen X, Qiu Z. The association between lactate dehydrogenase to serum albumin ratio and the 28-day mortality in patients with sepsis-associated acute kidney injury in intensive care: a retrospective cohort study. Ren Fail 2023;45(1):2212080. doi:10.1080/0886022X.2023.2212080, PMID:37194715.
- [68] Romanovsky A, Morgan C, Bagshaw SM. Pathophysiology and management of septic acute kidney injury. Pediatr Nephrol 2014;29(1):1–12. doi:10.1007/s00467-013-2427-6, PMID:23400860.
- [69] Atreya MR, Cvijanovich NZ, Fitzgerald JC, Weiss SL, Bigham MT, Jain PN, et al. Prognostic and predictive value of endothelial dysfunction biomarkers in sepsis-associated acute kidney injury: risk-stratified analysis from a prospective observational cohort of pediatric septic shock. Crit Care 2023;27(1):260. doi:10.1186/s13054-023-04554-y, PMID:37400882.